

**PATENT APPLICATION
5660-01102**

"EXPRESS MAIL" MAILING LABEL

NUMBER: 324271045
DATE OF DEPOSIT: July 11, 2003

I hereby certify that this paper is being deposited with the United States Postal Service "EXPRESS MAIL POST OFFICE TO ADDRESSEE" service under 37 C.F.R. 1.10 on the date indicated above and is addressed to: Commissioner for Patents, Alexandria, VA 22313-1450.



Derrick Brown

MODEL OF DENTAL CARIES

By:

Thomas D. Marshall

Attorney Docket No.: 5660-01102

Eric B. Meyertons
Meyertons, Hood, Kivlin, Kowert & Goetzl, P.C.
P.O. Box 398
Austin, Texas 78767-0398
Ph: (512) 853-8800

PRIORITY CLAIM

5 This application claims priority to U.S. Provisional Application No. 60/395,432
entitled “MODEL OF DENTAL CARIES” filed July 11, 2002. The above-referenced
provisional application is hereby incorporated by reference as if fully set forth herein.

BACKGROUND

10

1. Field of Invention

15 The present invention generally relates to a model of dental caries. Embodiments of the
invention relate to a model of dental caries that simulates dental decay in color, texture, and/or
tenacity.

2. Description of Related Art

20 Traditionally, pre-clinical dental students learn to treat various classifications of
decayed teeth by working on pristine resin teeth (teeth without defects) arranged in
simulated upper and lower resin jaws. In some cases, students may work with resin teeth
with defects. These defects may include various defects caused by human dental decay
(caries). Some resin teeth have defects but contain no simulated decay material. Other
resin teeth for pre-clinical use contain simulated decay material. Currently, simulated
25 carious resin teeth are produced by cutting away the outer layer of a resin crown, cutting
defects into the prepared tooth, inserting simulated carious material into the defects, and
then covering the prepared tooth with a resin crown. No opening is provided into the
dental cavity. Furthermore, the simulated carious material may be a heterogeneous
substance with unrealistic color, texture, and/or tenacity.

SUMMARY

A model of dental caries may include simulated decay material in a cavity in an
5 artificial tooth. In an embodiment, the artificial tooth may be made of resin, such as melamine resin. The cavity may include an opening that extends from an exterior surface of the tooth to an interior of the tooth. In some embodiments, an opening may extend at least to the dentinoenamel junction of the tooth. The opening may include at least one groove that extends from the opening along at least a portion of the dentinoenamel
10 junction. A groove may mimic progression of decay along the dentinoenamel junction. In an embodiment, one or more grooves are formed with a wheel bur. The cavity may be cleaned with an etching solution. The etching solution may be a liquid or a gel. In an embodiment, the etching solution may contain acid, such as phosphoric acid and/or orthophosphoric acid. The etching solution may be rinsed with water and the tooth may
15 be dried. A bonding substance, such as a liquid resin, may be applied to the walls of the cavity. In certain embodiments, the bonding substance may be chemically, light, and/or heat cured.

In various embodiments, a cavity may contain simulated decay material. The
20 simulated decay material may be in an opening and/or in at least one groove in the tooth. In an embodiment, the simulated decay material is a homogenous substance. The simulated decay material may consist of a curable resin and a porous substance. The curable resin may be chemically cured or cured with light and/or heat. In an embodiment, the porous substance may include pumice. In an embodiment, the
25 simulated decay material may include coloring, such as food coloring. The simulated decay material may simulate human dental decay in color, texture, and/or tenacity. The simulated decay material may be detectable with caries detecting stain. In some embodiments, the simulated decay material may be radio-opaque. In other embodiments, the simulated decay material may mimic the radiolucence of naturally occurring decay
30 material.

BRIEF DESCRIPTION OF THE DRAWINGS

Advantages of the present invention will become apparent to those skilled in the art with the benefit of the following detailed description and upon reference to the 5 accompanying drawings in which:

FIG. 1 depicts an opening in a resin tooth.

FIG. 2 depicts a wheel bur positioned to cut a groove in a resin tooth.

10

FIG. 3 depicts a schematic of a sagittal section of a cavity with an opening and a groove cut along a dentinoenamel junction of a resin tooth.

15

FIG. 4 depicts simulated decay material in a groove along the dentinoenamel junction.

FIG. 5 depicts a schematic of a sagittal section of simulated decay material in a groove cut along a dentinoenamel junction of a resin tooth.

20

FIG. 6 depicts simulated decay on a side of an opening in a resin tooth.

FIG. 7 depicts simulated decay material in an opening in a resin tooth.

25

FIG. 8 depicts a schematic of a sagittal section of simulated decay material in a cavity in a resin tooth.

FIG. 9 depicts simulated occlusal decay.

While the invention may be susceptible to various modifications and alternative 30 forms, specific embodiments thereof are shown by way of example in the drawings and will herein be described in detail. The drawings may not be to scale. It should be

understood, however, that the drawings and detailed description thereto are not intended to limit the invention to the particular form disclosed, but to the contrary, the intention is to cover all modifications, equivalents, and alternatives falling within the spirit and scope of the present invention as defined by the appended claims.

5

DETAILED DESCRIPTION

FIG. 1 depicts artificial tooth 100 with opening 102. In an embodiment, tooth 100 may be made of resin, such as melamine resin, or any other material approved by the
10 American Dental Association for use as artificial teeth for use in dental education and for use in state and/or regional board examinations. Opening 102 may extend from the enamel or exterior surface to an interior of the tooth. In some embodiments, opening 102 may extend from the enamel or exterior surface to at least an artificial dentinoenamel junction. In certain embodiments, opening 102 extends beyond the artificial
15 dentinoenamel junction into artificial dentin.

In FIG. 2, wheel bur 104 is shown inserted into opening 102 of tooth 100. In an embodiment, wheel bur 104 may be used to cut one or more grooves extending from the opening along the dentinoenamel junction. In other embodiments, other instruments may
20 be used to cut one or more grooves along the dentinoenamel junction. As used herein, a “groove” is an opening including, but not limited to, a valley, a notch, etc. In some embodiments, a groove may be formed around a full circumference of the opening. In other embodiments, one or more grooves may be formed around a portion or portions of the circumference of the opening. The grooves are cut to simulate progression of human
25 dental decay as the decay spreads from the opening in the enamel along the dentinoenamel junction. In certain embodiments, grooves may be cut in various locations in the artificial tooth. One or more grooves may be completely or partially filled with simulated decay material to form various embodiments of a tooth with caries. A groove with simulated decay material may provide a representative model of caries along the
30 dentinoenamel junction of a tooth.

A schematic of a sagittal section of tooth 100 is shown in FIG. 3. “Sagittal section” is used herein to refer to a longitudinal cross section. Artificial enamel layer 106 and artificial dentin 108 form an interface at dentinoenamel junction 110. Groove 112 is
5 shown extending from opening 102 along dentionenamel junction 110. In an embodiment, a cavity in tooth 100 may include opening 102. Opening 102 may include one or more grooves 112. Grooves 112 simulate the spread of dental decay along dentionenamel junction 110. In the embodiment of FIG. 3, opening 102 extends beyond dentinoenamel juction 110 into dentin 108. FIG. 4 depicts simulated decay material 114
10 substantially in groove 112 extending from opening 102 along the dentinoenamel junction of tooth 100.

FIG. 5 depicts a schematic of a sagittal section of tooth 100 with a cavity including opening 102 and groove 112. In the embodiment of FIG. 5, opening 102
15 extends through artificial enamel 106, beyond dentinoenamel junction 110, and into dentin 108. At least a portion of groove 112 may contain simulated decay material 114.

FIG. 6 depicts an embodiment of simulated decay material 114 in opening 102 of tooth 100. In the embodiment of FIG. 6, simulated decay material 114 is located
20 substantially along a side of opening 102 in tooth 100. FIG. 7 depicts another embodiment of simulated decay material 114 in an opening in tooth 100. In the embodiment of FIG. 7, simulated decay material 114 extends beyond the dentinoenamel junction and substantially fills the opening in tooth 100. FIG. 8 depicts a schematic of a sagittal section of simulated decay material 114 in opening 102 in tooth 100. In the
25 embodiment shown in FIG. 8, simulated decay material 114 fills part of opening 102 in tooth 100. In this embodiment, at least part of groove 112 along dentinoenamel junction 110 between artificial enamel 106 and dentin 108 contains simulated decay material 114. FIG. 9 depicts an occlusal view of tooth 100 with simulated occlusal decay shown by simulated decay material 114.

FIGS. 4-9 show various embodiments of defects and simulated decay in artificial teeth. Other embodiments may include various types and sizes of defects, as well as various extents of decay. In some embodiments, decay in a model tooth may be limited to a small portion of an opening in the tooth. In other embodiments, a model tooth may
5 be made to exhibit extensive decay in a pulpal direction and/or along grooves in the dentinoenamel junction. Openings may be formed on any exterior surface of an artificial tooth. The composition of simulated decay material may be chosen to simulate the color, texture, and/or tenacity of dental caries at various stages of decay.

10 In an embodiment, an artificial tooth with simulated decay may be formed from a pristine resin tooth. A defect of a desired character may be formed in the tooth with an appropriate cutting or drilling instrument. A groove may be formed in the tooth along the dentinoenamel junction, as shown in FIG. 2. An etching solution, such as a gel or a liquid, may be applied to surface defined by a cavity, including a surface defined by an
15 opening and/or a surface defined by one or more grooves. The etching solution may contain an acid. In certain embodiments, the etching solution may include phosphoric acid and/or orthophosphoric acid. Specifically, the etching solution may be a solution of orthophosphoric acid (e.g., a solution containing 20-50% acid in water or a solution containing 37% acid in water). The etching solution may be applied to a surface defined
20 by the cavity with a brush or similar applicator with a felt tip or sponge tip. The etching solution may be allowed to remain on the resin tooth for a period of at least 5 seconds to 60 seconds or longer, such as 10-30 seconds or 15-20 seconds, depending on the composition and concentration of the etching solution. The etching solution may be rinsed from the tooth (with, e.g., water). The tooth may be dried.
25

A bonding substance, such as a liquid resin, including unfilled or lightly filled resins, may be applied to a surface defined by the cavity. The bonding substance may wet the tooth. In an embodiment, the bonding substance may adhere to the resin tooth and bond to synthetic resin. In various embodiments, the bonding substance may be a
30 commercial dentin bonding agent or a mixture of methyl methacrylate and methylene chloride. In other embodiments, the bonding substance may be any substance capable of

wetting and adhering to resin teeth. The bonding substance may be applied with a brush or similar applicator with a felt tip or sponge tip. Preparation of the tooth for application of simulated decay material may include curing the bonding substance with heat or light (e.g., visible, ultraviolet, or infrared light). In an embodiment, the bonding substance
5 may be cured with visible light with a wavelength of about 500 nm and an output of about 300 milliwatts/square centimeter for about 10 seconds to about 30 seconds, or about 20 seconds. A light source such as the ESPE-Highlight manufactured by 3M (St. Paul, MN) or other similar light or lights may be used.

10 Simulated decay material may be applied to a cavity, including the opening and/or one or more grooves. Some naturally occurring dental caries may have a soft outer layer and increasingly harder inner layers. Other naturally occurring dental caries may have a substantially uniform hardness. Applying simulated decay material of a uniform composition may result in simulated dental caries of uniform hardness. Layering of
15 simulated decay material of various compositions may result in simulated dental caries with layers of varying hardness. In some embodiments, an inner layer of simulated decay material may have a greater hardness than an outer layer of simulated decay material. In other embodiments, an outer layer of simulated decay material may have a greater hardness than an inner layer of simulated decay material.

20 Simulated decay material may include a resin. In an embodiment, the resin is non-hardening. The resin may be any material approved for use by the American Dental Association of suitable texture (i.e., consistency) to simulate human dental caries. The resin may be heat cured and/or chemically cured. In an embodiment, the resin may be
25 cured with light. The resin may be radio-opaque or radiotranslucent. Tempit® L/C, manufactured by Centrix Incorporated (Shelton, Connecticut), a resin which tends to stay pliable after curing, may be used in simulated decay material. Properties of the resin may allow the simulated decay material to be removed from a tooth with an excavator or other dental instrument. In an embodiment, caries detecting stain may be used to check for
30 complete removal of simulated decay material in a pre-clinical setting.

Simulated decay material may include a porous substance. The porous substance may include pumice or similar material. The porous substance may be approved for use by the American Dental Association. Nupro® T Prophylaxis Paste manufactured by Dentsply Preventive Care (York, Pennsylvania) may be used as the porous substance. In some embodiments, simulated decay material may contain substantially equal portions of resin and porous substance. In other embodiments, depending on the desired texture of the simulated decay material, the ratio of resin to porous substance may range from about 1:5 parts by volume resin to porous substance to about 5:1 parts by volume resin to porous substance. A 1:1 volume ratio of resin to porous substance may result in simulated decay material with an intermediate texture (e.g., intermediate hardness). A greater volume of resin than porous substance may result in simulated decay material with a harder texture. A greater volume of porous substance than resin may result in simulated decay material with a softer texture. In some embodiments, a homogeneous simulated decay material may be used in to form artificial caries of uniform hardness. In other embodiments, homogeneous simulated decay materials of various compositions may be layered in cavities of artificial teeth to create simulated caries of varying hardness.

Simulated decay material may be detected by a caries detecting stain, such as Cari-D-Tect, manufactured by Gresco Products, Inc. (Stafford, Texas). In some embodiments, simulated decay material may be radio-opaque. In other embodiments, simulated decay material may be radiotranslucent. Upon X-ray, radiotranslucent simulated decay material may appear substantially the same as naturally occurring decay material.

In an embodiment, coloring may be added to the resin/porous substance mixture to simulate the color of human dental caries. The coloring may be food coloring or any coloring approved for use by the American Dental Association or the Federal Drug Administration. Coloring may be added to achieve a desired coloring of simulated decay

material. In various embodiments, the desired color may range from dentin color to tan, to shades of brown, gray, and/or black.

In some embodiments, the simulated decay material may be chemically cured or
5 heat cured. In certain embodiments, the simulated decay material may be cured with light, such as visible light with a wavelength of about 500 nm, to achieve the desired texture and tenacity (i.e., degree of bonding to the tooth). In an embodiment, the simulated decay material may be cured with visible light with a wavelength of about 500 nm and an output of about 300 milliwatts/square centimeter for about 10 seconds to about
10 70 seconds, about 20 seconds to about 60 seconds, about 30 seconds to about 50 seconds, or about 40 seconds. Simulated decay material of a darker color may require a longer curing time than simulated decay material of a lighter color.

Simulated decay material in a cavity of an artificial tooth may be cured with one
15 or more applications of heat and/or light. In some embodiments, simulated dental caries may be formed by curing simulated decay material of uniform hardness with one application of heat and/or light. In certain embodiments, simulated dental caries may be formed by successively curing layers of simulated decay material. Successively curing layers of simulated decay material of uniform hardness may result in simulated dental
20 caries of uniform hardness. Successively curing layers of simulated decay material of varying hardness may result in simulated dental caries with layers of varying hardness.

The model of dental caries described in the above embodiments is a realistic representation of human dental decay. Advantages of certain embodiments of this model
25 include a substantially homogeneous simulated decay material that mimics human carious dentin in color, texture, and tenacity. Another advantage of certain embodiments is the sensitivity of the simulated decay material to caries detecting stain. Another advantage of the model of dental caries is the presence of one or more grooves along the dentinoenamel junction, replicating the progression of dental caries in the human
30 condition. Another advantage of certain embodiments of the model of dental caries

described above is the low cost and ready availability of materials needed to incorporate simulated decay into artificial teeth.

Further modifications and alternative embodiments of various aspects of the
5 invention will be apparent to those skilled in the art in view of this description.

Accordingly, this description is to be construed as illustrative only and is for the purpose
of teaching those skilled in the art the general manner of carrying out the invention. It is
to be understood that the forms of the invention shown and described herein are to be
taken as examples of embodiments. Elements and materials may be substituted for those
10 illustrated and described herein, parts and processes may be reversed, and certain features
of the invention may be utilized independently, all as would be apparent to one skilled in
the art after having the benefit of this description of the invention. Changes may be made
in the elements described herein without departing from the spirit and scope of the
invention as described in the following claims.

15